

How to take a picture when ordering an exoskeleton glove

In order to secure you the best fit of your exoskeleton glove, it's important that the following steps are done as precisely as possible.

Step #1 Close your hand softly. Use a **black** white-board marker and with a circular motion make dots of 3-5 mm in diameter (as seen in figure 3) on each finger, at the positions seen in figure 1.




Figure 1 - Where to mark the center of your joints on each finger

Step #2 Place your hand in the middle of a white piece of A4 paper, on a table which is NOT white

Step #3 Position your phone vertically above the center of the hand (See figure 2). If you stand up with the piece of paper on a desk, the height of your phone should be between your elbow and shoulder. Make sure all corners of the piece of paper are visible!




Figure 2 - Wrong and correct position of the camera relative to the background

Step #4 Spread your fingers just enough to see the background clearly, as illustrated in figure 3.




Figure 3 - Amount of finger spread: Wrong on the left and Correct on the right

Step #5 take the picture with the FLASH ON and with a steady hand.